

A decorative graphic on the left side of the slide, featuring a vertical column of seven hexagons. From top to bottom: a teal hexagon, a grey hexagon, a green hexagon with a microscope icon, a grey hexagon with a DNA double helix icon, a teal hexagon with a document icon, a green hexagon with a medicine bottle icon, and a teal hexagon. The background of the slide is a light grey honeycomb pattern.

Comprehensive Genomic Profiling with FMI – bringing value to patients

Vasanti Natarajan, PhD

Medical lead Molecular Information

Region Europe

25 th October 2016



Why should physicians consider profiling?



CGP from Foundation Medicine – *understanding the difference*




Foundation Medicine's experience – *translating into clinical benefit*

Cancer is a disease of the genome


Our understanding of cancer has been evolving




Etiology




Genomic Alteration



Malignancy



Many genetically driven characteristics— *many therapeutic options*




Lung Adenocarcinoma:

Moving from one disease to multiple disease types by molecular alterations that require distinct tx plans



2000





2015

2015 NSCLC NCCN guidelines recommend broad molecular profiling for the following biomarker/drug associations:

Biomarker	Drugs
EGFR mutations	erlotinib, gefitinib, afatinib
ALK rearrangements	crizotinib
BRAF V600E	vemurafenib*, dabrafenib*
MET amplifications	crizotinib
ROS1 rearrangements	crizotinib
HER2 mutations	trastuzumab*, afatinib
RET rearrangements	cabozantinib*

* Drugs not approved for lung cancer




Treating Lung Cancer patients based on their tumour profiling results improves outcomes



Meta-analysis of Phase II studies – 32 149 patients

Meta-analysis of 570 Phase II, single-agent studies (including total of 32,149 patients) studying the impact of personalized and targeted treatment strategies in diverse cancer types

RR, PFS and OS from pooled meta-analyses



Arms included (pooled and meta-analyses):
RR: Personalized n = 112 not personalized n = 526
PFS: Personalized n = 86 Not personalized n = 444
OS: Personalized n = 49 Not personalized n = 392

Personalized treatment strategies, across malignancies, were **independent predictors of better outcomes and fewer deaths** from treatment toxicity than non-personalized therapies

“Matched therapy using genomic markers offers better outcomes than using protein biomarkers”

“Matched therapies are associated with better outcomes than non-matched therapies”

Schwaederle, M., et al. (2015), J Clin Oncol 33(32):3817-25.

OS: Overall survival; PFS: Progression-free survival; RR: Response rate; CI: Confidence interval

Profiling guidelines

NCCN Guidelines® now recommend “broad molecular profiling” for advanced NSCLC patients

Genomic Alterations (i.e. driver event)	Available targeted agents with activity against driver event in lung cancer
EGFR mutations	erlotinib, gefitinib, afatinib
ALK rearrangements	crizotinib, ceritinib
HER2 mutations	trastuzumab, afatinib
BRAF V600E mutations	vemurafenib, dabrafenib
MET amplification	crizotinib
ROS1 rearrangements	crizotinib
RET rearrangements	cabozantinib

NCCN Guidelines® Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines®) for Non-Small Cell Lung Cancer V.4.2016. © National Comprehensive Cancer Network, Inc. 2016. All rights reserved. Accessed March 13, 2016. To view the most recent and complete version of the guideline, go online to NCCN.org. NATIONAL COMPREHENSIVE CANCER NETWORK®, NCCN®, NCCN GUIDELINES®, and all other NCCN Content are trademarks owned by the National Comprehensive Cancer Network, Inc. now recommends “broad molecular profiling” for advanced NSCLC patients

Profiling initiatives

Investigating the potential to match treatments to genomic alterations across tumor types



CAPTUR Canadian Profiling and targeted Utilization trial



DRUP The Drug Rediscovery Protocol



TAPUR Targeted Agent and Profiling Utilization Registry Study



And more...



Initiatives to decipher which patients respond to which therapies, irrespective of in which tumor type the therapies are approved in

http://www.nature.com/nm/journal/v22/n5/fig_tab/nm.4089_T1.html



Why should physicians consider profiling?



Comprehensive genomic profiling from Foundation Medicine



How to benefit from Foundation Medicine's experience

Foundation Medicine

Pioneer and leader in molecular information

- Founded 2010 in Cambridge, MA, USA
- Proprietary molecular information platform
- First to market comprehensive genomic profiling solutions for cancer
- 90,000+ clinical cases profiled
- 30+ pharmaceutical clinical trial partners
- Roche collaboration for R&D and commercialization outside USA



Foundation Medicine offers two solutions

FoundationOne® and FoundationOne® Heme



Comprehensive: Detect all classes of genomic alterations



FOUNDATIONONE™

Coding regions of 315 genes

Introns of 28 genes

Known as drivers of solid tumors

Frampton G, et al. Nature Biotech, 2013, 31, 1023-34



FOUNDATIONONE™
Heme

DNA sequences of 405 genes

RNA sequences (cDNA) of 265 genes

Commonly altered in hematologic malignancies (leukemia, lymphoma and myeloma) and

Jie H et al. Blood 2016


How does FoundationOne work?



A process that follows standard operating processes

Pre-Analytic Process (Pre-Sequencing)

Post-Analytic Process (Post-Sequencing)



Powered by 20+ bioinformaticians and genomic scientists who optimize state-of-the-art algorithms to report the most clinically relevant information for a patient



Why should physicians consider profiling?



CGP from Foundation Medicine – *understanding the difference*



How to benefit from Foundation Medicine's experience

Types of genomic alterations driving tumor growth

Limitations of traditional and hotspot testing



BRAF V600E
(BRAF inhibitor)
Base Substitutions



EGFR Exon 19 deletion
(EGFR inhibitor)
Insertions and Deletions



HER2 amplification
(HER2 inhibitor)
Copy Number Alterations



ALK fusion
(ALK inhibitor)
Rearrangements

Test	Detects	Can Miss
IHC	Protein expression	Any alteration not known of ahead of time
FISH	Copy number alterations, Rearrangements	Insertions & deletions, Substitutions
Hot Spot NGS*	Substitutions	Insertions & deletions, Copy number alterations, Rearrangements

*Suh J et al. (2016) *Oncologist*. <http://dx.doi.org/10.1634/theoncologist.2016-0030>.


Foundation Medicine finds more targets

Completely sequencing genes enables detection of novel alterations missed by hot spot testing

Example: EGFR gene

Hot spot tests detect selective alterations in selective parts of the EGFR gene*

....When there is an insertion/deletion or rearrangement that removes one of the primer sites, hot spot tests will not amplify the region or detect the alteration



FoundationOne detects all genomic alterations across the entire EGFR gene

= Mutations not detected by hot spot

* Meric-Bernstam F et al. (2015) J Clin Oncol 33:2753-2762

Foundation Medicine finds more targets

While hot spot tests can miss alterations...



MULTI-GENE "hot spot" TEST




Foundation Medicine finds more targets


...Comprehensive Genomic Profiling identifies all four classes of alterations with validated performance

COMPREHENSIVE GENOMIC PROFILING

MET exon 14 splice
ALK inhibitor



EGFR L858R
EGFR inhibitor




ROS1 fusion
ALK inhibitor



MET amplification
ALK inhibitor



EGFR exon 19 deletion
EGFR inhibitor



Validated performance
published in peer-
reviewed journal*

Base substitutions:
Sensitivity: >99%
PPV: >99%

Rearrangements:
Sensitivity: ≥90%
PPV: >99%

**Copy number
alterations:**
Sensitivity: >95%
PPV: >99%

Insertions/deletions:
Sensitivity: >97%
PPV: >99%

* Frampton G et al. (2013) Nature Biotech 31, 1023-34



Why should physicians consider profiling?



Comprehensive genomic profiling from Foundation Medicine



How to benefit from Foundation Medicine's experience




Improve profiling of NSCLC patients

FoundationOne finds more alterations associated with NCCN guidelines than single gene or hot spot NGS

SINGLE GENE TESTING

misses up to 35% of ALK rearrangements by FISH¹ and 17% of EGFR alterations by hot spot test²



Hot spot NGS

Up to 50% of targetable alterations can be missed without supplemental FISH³



detects all four classes of NSCLC clinically relevant alterations³ and genetic biomarkers⁴ included in the NCCN Guidelines[®]



1. Ali S et al. (2016) *Oncologist*. doi:10.1634/theoncologist.2015-0497.
2. Schrock AB et al. (2016) *Clin Cancer Res*. Mar 1. pii: clincanres.1668.2015.
3. Suh J et al. (2016) *Oncologist*. <http://dx.doi.org/10.1634/theoncologist.2016-0030>
4. NCCN Guidelines[®] Referenced with permission from the NCCN Clinical Practice Guidelines in Oncology (NCCN Guidelines[®]) for Non-Small Cell Lung



Clinical utility of finding more alterations with FoundationOne

NSCLC patients can benefit from targeted therapies



17%

of EGFR exon 19 deletions missed by hotspot tests¹

35%

of ALK-rearranged cases missed by FISH³

75%

of NSCLC patients with EGFR exon 19 deletions can respond to EGFR tyrosine kinase inhibitors, with median OS > 1 year²

80%

of ALK-rearranged patients identified by FoundationOne respond to ALK inhibitor crizotinib³

Patient case: EGFR/ALK negative

Identification of complex fusion led to treatment/response


Patient Information



- 43-year-old male
- Never-smoker

At Presentation

- Pericardial tamponade
- No detection of *EGFR* mutation; atypical FISH staining for *ALK*



Atypical pattern of double 3'ALK signals (red) fused with 5'ALK signal (green)

Diagnosis


- **Metastatic NSCLC** with a pericardial tamponade

Treatment status

- Disease **progression** despite 4 cycles of cisplatin/pemetrexed

FoundationOne® analysis and subsequent treatment

- Identification of complex EML4-ALK fusion separated by genomic shards
- Initiation of treatment with crizotinib



Chest PET-CT before (*) and after 4 months (**) crizotinib treatment showing a significant shrinkage of the primary lesion

Patient had a rapid response to crizotinib treatment; **75% shrinkage** of primary lesion (RECIST) after 4 months of treatment

FoundationOne may lead to improved outcomes

Studies show potential in other tumor types, ability to impact physician decisions



Triple-negative breast cancers

Targeted therapy based on tumor genomic alterations

Improved response (33% vs. 8%, $p=0.018$) & **longer progression-free survival**¹ (6.4 vs. 1.9 months; $p=0.001$)



Rare/refractory gynecological cancers

Targeted therapy based on tumor genomic alterations

Radiologic response or stability in 64% of patients²



Breast cancer

Targeted therapy based on tumor genomic alterations

41% of treatment decisions influenced by FoundationOne³



Targeting *ERBB2* Mutations in Metastatic Breast Cancer

Human Cancer Biology

Clinical
Cancer
Research

Relapsed Classic E-Cadherin (*CDH1*)–Mutated Invasive Lobular Breast Cancer Shows a High Frequency of *HER2* (*ERBB2*) Gene Mutations

Jeffrey S. Ross^{1,2}, Kai Wang², Christine E. Sheehan¹, Ann B. Boguniewicz¹, Geoff Otto², Sean R. Downing², James Sun², Jie He², John A. Curran², Siraj Ali², Roman Yelensky², Doron Lipson², Gary Palmer², Vincent A. Miller², and Philip J. Stephens²

- Responses reported with both antibody therapeutics and kinase inhibitors
- 38% response rate in *ERBB2* mutated BC to kinase inhibitor at SABCS
- High frequency (> 30%) of *ERBB2* mutations in *CDH1* mutated relapsed ILC

RESEARCH ARTICLE

Activating *HER2* Mutations in *HER2* Gene Amplification Negative Breast Cancer

Ron Bose^{1,2}, Shyam M. Kavuri¹, Adam C. Searleman¹, Wei Shen¹, Dong Shen³, Daniel C. Koboldt³, John Monsey¹, Nicholas Goel¹, Adam B. Aronson¹, Shunqiang Li^{1,2}, Cynthia X. Ma^{1,2}, Li Ding^{1,2,3,4}, Elaine R. Mardis^{2,3,4}, and Matthew J. Ellis^{1,2}

JNCCN

The official journal of the

NCCN


Nat
Cor
Car
Net

HOME | CURRENT ISSUE | PAST ISSUES | SUPPLEMENTS | CME | CONTACT US | SUBSCRIPTIONS

HER2-Mutated Breast Cancer Responds to Treatment With Single-Agent Neratinib, a Second-Generation *HER2*/*EGFR* Tyrosine Kinase Inhibitor

Noa Efrat Ben-Baruch, MD^a, Ron Bose, MD, PhD^{b,c}, Shyam M. Kavuri, PhD^{b,c,d}, Cynthia X. Ma, MD, PhD^{b,c} and Matthew J. Ellis, MB, BChir, PhD, FRCP^{b,c,d}

Response of a HER2 FISH/IHC Negative Cutaneous Adnexal Carcinoma with an *ERBB2* S310f Mutation to anti-HER2 Targeted Therapy



Vornikova O et al. *The Oncologist* 2014;19:1006-1007

Why consider profiling with Foundation Medicine?

- Profiling has been shown to **improve outcomes for patients** with lung cancer or considering clinical trials, while evidence is evolving in additional indications¹⁻³
- Foundation Medicine's profiling services are designed to **capture all four types of genomic alterations** which single gene and hotspot NGS testing can miss
- Proprietary **bioinformatics** have been optimized over 90,000 cases to call alterations
- These alterations are delivered in a **comprehensive report** which describes potential therapies, trials, and the latest clinical literature to inform physician's decisions
- Evidence has shown **FoundationOne detects alterations in patients that are pan-negative with single gene panels**⁴⁻⁶, and in some indications can improve outcomes^{8-4, 7-8}



Profiling with FoundationOne finds more clinically-relevant alterations and can lead to better patient outcomes

Doing now what patients need next